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## 12 Ihre Notizen

## Table of contents

1	introduction .....	49
2	certification of <i>piezosystem jena</i> .....	49
3	declaration of conformity .....	49
4	purchased part package .....	50
5	instructions for using piezoelectric elements and power supplies .....	50
6	safety instructions .....	51
6.1	installation, power supply .....	52
6.2	operation.....	52
6.3	maintenance and inspection .....	53
6.4	environmental conditions .....	53
7	instructions for checking the function of the system / quick start .....	53
8	how to operate the digital system d-Drive®.....	54
8.1	common introduction .....	54
8.2	system housing 84TE with power supply .....	55
8.2.1	function .....	55
8.2.2	technical data system housings .....	55
8.3	power amplifier EVD50, EVD125, EVD300 .....	56
8.3.1	function .....	56
8.3.2	MOD/MON .....	59
8.3.3	pinning .....	61
8.3.4	technical data EVDxxx .....	62
8.4	interface EDS1, EDS2 .....	63
8.4.1	function EDS2 .....	63
8.4.2	menus EDS2 .....	64
8.4.3	pinning EDSx.....	65
8.4.4	technical data EDSx.....	66
8.5	communication & commands.....	66
8.5.1	communication via RS232 .....	66
8.5.2	driver installation of USB interface.....	67
8.5.3	commands .....	73
8.5.4	status register.....	77
8.6	function generator .....	77
8.7	output of trigger signals.....	78
8.8	scan function .....	81
8.9	data recorder .....	83
9	handling .....	85
10	controller adjustment.....	86
11	troubleshooting.....	87
11.1	error register .....	88
12	your notes.....	89

## 1 introduction

This manual describes the **d-Drive**<sup>®</sup> piezo digital system from **piezosystem jena**. You will also find additional information regarding piezoelectric products.

Definition: All systems from **piezosystem jena**, such as electronics, actuators and optical systems, are called “units”.

If you have any problems please contact the manufacturer of the system: **piezosystem jena**, Stockholmer Strasse 12, 07747 Jena. Phone: +49 36 41 66 88-0

## 2 certification of **piezosystem jena**

The company, piezosystem jena GmbH, has worked according to a DIN EN ISO 9001 certified quality management system since 1999. Its effectiveness is verified and proven by periodic audits by the TÜV.

This instruction manual includes important information about using piezo actuators. Please take the time to read this information. Piezo positioning systems are mechanical systems that offer the highest precision. Correct handling guarantees the maintenance of this precision over a long period of time.

## 3 declaration of conformity

The CE Declaration of Conformity is available on our website [www.piezosystem.com](http://www.piezosystem.com).

#### 4 purchased part package

Please check the completeness of the delivery after receiving the shipment:

- **d-Drive**<sup>®</sup> digital system
- power cord
- RS232 cable
- USB cable
- instruction manual
- CD-ROM with driver, software, and instruction manual

#### 5 instructions for using piezoelectric elements and power supplies

- Piezoelectric actuators from **piezosystem jena** are controlled by voltages of up to 150 V. These values can be quite hazardous. Therefore, read the installation instructions carefully and ensure that only authorized personnel handle the power supply.
- After transportation, piezoelectric actuators should be allowed to adapt to room temperature for approximately 2 hours before being switched on.
- Piezoelectric actuators are made from ceramic materials with and without metallic casings. The piezo-ceramic is a relatively brittle material. This should be noted when handling piezoelectric actuators. All piezo elements are sensitive to bending or shock forces. Due to the piezoelectric effect, piezo actuators can generate electrical charges by changing the mechanical load or the temperature, or by such actions as those described above.
- Piezoelectric actuators are able to work under high compressive forces. Only actuators with a pre-load can be used under tensile loads (these tensile forces must be less than the pre-load, given in the data sheet). Please note that the acceleration of the ceramic material (e.g., caused by fall down, discharging, or high dynamic application) will occur.
- After excitation of the actuators by a voltage in the upper control range, the ceramic will move and generate an opposite high voltage after disconnection.
- Heating of the ceramic material will occur during dynamic operation and is caused by structure conditional loss processes. This may cause failure if the temperature exceeds the specified values cited below. With increasing temperature, up to the Curie temperature  $T_C$  (usual values approx. 140°C - 250°C), the piezoelectric effect disappears. We recommend working in temperatures up to  $T_C/2$  (normally up to 80°C).
- Piezoelectric actuators, such as stacks or other devices, work electrically as a capacitor. These elements are able to store electrical energy over a long period of time (up to some days) and the stored energy may be dangerous.
- If the actuator remains connected to the drive electronics, it will be unloaded within a second after shutdown and quickly reaches harmless voltage values.
- Piezo actuators can only generate voltages by warming or cooling (caused by the longitudinal change). The discharge potential should not be ignored due to the inner capacitance. This effect is insignificant at room temperature.
- Piezo actuators from **piezosystem jena** are adjusted and glued. Any opening of the unit will cause misalignment or possible malfunction, and will result in the invalidation of the guarantee.
- Please only use original parts from **piezosystem jena**.
- Please contact **piezosystem jena**, or your local representative, if there are any problems with your actuator or power supply.

Caution!

Shock forces may damage the built-in ceramic elements. Please avoid such forces, and handle the units with care, otherwise the guarantee will become invalid.

## 6 safety instructions

### Icons:

**RISK OF ELECTRIC SHOCK!** Indicates that a risk of electric shock is present and the associated warning should be observed.

**CAUTION! REFER TO OPERATOR'S MANUAL** – Refer to your operator's manual for additional information, such as important operating and maintenance instructions.

### **RISK OF ELECTRIC SHOCK!**

- Do not open the units! There are no user serviceable parts inside and opening or removing covers may expose you to dangerous shock hazards or other risks. Refer all servicing to qualified service personnel.
- Do not spill any liquids into the cabinet or use the units near water.

### **CAUTION!**

- Allow adequate ventilation around the units so that heat can properly dissipate. Do not block ventilated openings or place the units near a radiator, oven, or other heat sources. Do not put anything on top of the units except those that are designed for that purpose (e.g. actuators).
- Only work with the units in a clean and dry environment! Only specially prepared units (e.g. actuators) can work under other conditions!
- Please only use original parts from **piezosystem jena**. **piezosystem jena** does not give any warranty for damages or malfunction caused by additional parts not supplied by **piezosystem jena**. Additional cables or connectors will change the calibration and other specified data. This can change the specified properties of the units and cause them to malfunction.
- Piezo elements are sensitive systems capable of the highest positioning accuracy. They will only demonstrate their excellent properties if they are handled correctly! Please mount them properly only at the special mounting points.

Immediately unplug your unit from the wall outlet and refer servicing to qualified service personnel under the following conditions:

- when the power supply cord or plug is damaged
- if liquid has been spilled or objects have fallen into the unit
- if the unit has been exposed to rain or water
- if the unit has been dropped or the housing is damaged

## 6.1 installation, power supply

### RISK OF ELECTRIC SHOCK

- Do not insert or unplug the power plug with wet hands, as this may result in electrical shock.
- Do not install in rooms where inflammable substances are stored. If flammable substances come into contact with electrical parts inside, it could result in fire or electrical shock.
- Do not damage or modify the power cord. Also, do not place heavy objects on the power cord, or pull on or excessively bend it, as this could cause electrical damage and result in a fire or electrical shock.
- Always grasp the plug portion when unplugging the power cord. Pulling on the power cord may expose or snap the core wire, or otherwise damage the power cord. If the cord is damaged, this could cause an electricity leak and result in a fire or electrical shock.

### CAUTION!

- Do not use accessories other than the ones provided (e.g. power cord). Only plug the power cord into grounded power equipment and sockets.
- Do not place heavy objects on any cables (e.g. power cords, sensor cables, actuator cables, optical cables).
- Do not block ventilated openings or place the units near a radiator, oven, or other heat sources.
- Plug in the power cord completely so that it cannot loosen inadvertently.
- Leave sufficient space around the power plug so that it can be unplugged easily. If objects are placed around the power plug, you will be unable to unplug it in an emergency.
- Install the system so that the on/off-switch is easily accessible at all times.
- The power plug is the cut-off point to the main power supply.

## 6.2 operation

### RISK OF ELECTRIC SHOCK!

- Do not open the units! There are no user serviceable parts inside and opening or removing covers may expose you to dangerous shock hazards or other risks. Refer all servicing to qualified service personnel.
- Do not spill inflammable substances inside the voltage amplifier. If these items come into contact with an electrical component inside the voltage amplifier, it may result in a fire or electrical shock.

### CAUTION!

- If the voltage amplifier emits smoke, high heat, or unusual smells, immediately turn off the power switch and unplug the power plug from the outlet. Then contact our technical service.

### 6.3 maintenance and inspection

#### CAUTION!

- Before cleaning the exterior box of the voltage amplifier, turn off the power switch and unplug the power plug. Failure to do so may result in a fire or electrical shock.
- Clean the exterior box using a damp cloth that has been firmly wrung-out. Do not use alcohols, benzene, paint thinner, or other inflammable substances. If flammable substances come into contact with an electrical component inside the voltage amplifier, it may result in a fire or electrical shock.

### 6.4 environmental conditions

The amplifier can be used:

- indoors only
- at an altitude up to 2000 m
- at a temperature of 5 ... 35 °C
- at a relative humidity of 5 ... 95% (non-condensing)

The recommended environmental conditions:

- indoors only
- at an altitude up to 2000 m
- at a temperature of 20...22 °C
- at a relative humidity of 5 ... 80% (non-condensing)

## 7 instructions for checking the function of the system / quick start

When you open the package, please check to make sure all the necessary parts are included (see packing list) and nothing is damaged.

Check the electronics and the actuator for any visible damage:

- The top and bottom plates of the actuator (if it does not have another shape) should be parallel each to each other, and not contain any scratches.
- If there is any damage to the system please contact our local representative immediately!
- If the packaging material is damaged please confirm this with the shipping company.
- Ensure that the main voltage supplied in your country is the same as the one installed for the system (Check the voltage label on the backside of the power supply).
- The power switch should be in the off position.
- Connect the power supply.
- Connect the piezo actuator by using the D-SUB 15pin connector. Be sure that the cables are connected properly to the electronics.
- Switch on the power supply on the back side (Standby).
- Switch on the system by pressing the POWER button on the EDS. The green "READY" LED indicates the power supply is working.
- It is followed by a short self test of all parts.
- All of the yellow LEDs included in the EVD50 will blink. The green "ON" LED will shine permanently. If there is an actuator plugged into the D-SUB 15pin and the amplifiers are working, the LED's will switch off after about 3 seconds.

- The interface EDS2 finishes its self test by showing “main menu”, and then the channels of amplifiers present are displayed. If there are amplifiers in slot 1 and 2, the display in rows 1 and 2 shows the actuator name “X” (“Y”...), “ol” (open loop) and the actuator voltage “-20.000 V”. The rows 3...6 are blank because there are no amplifiers or actuators connected. If the system includes the EDS1 without display, you can check the status with the command “stat *Enter*” by using a terminal program via RS232.
- Turn the potentiometer to the right until the display shows about 130 V. The actuator is in the maximum position.
- Turn the potentiometer to the left till the display shows about -20 V. The actuator is in the minimum position.
- Switch on closed loop by pressing on the potentiometer. The “OL/CL” LED will shine yellow. On the EDS2 the display will change from “op” to “cl” and “-20.000 V” to “0  $\mu$ m”. Turn the potentiometer to the right until the maximum position is reached and the actuator does its maximum closed loop movement.
- Turn the potentiometer to the left until the minimum position (0  $\mu$ m) is reached and switch off the closed loop.
- Do this procedure with all connected actuators.
- Switch off the power supply by pressing the POWER button on the EDS2 longer than 2 seconds.
- If the system is switched off a long time (longer then 1 day), please switch off the switch on the back side of the system.
- The function check is complete.
- If there were errors during the procedure, read the *troubleshooting* chapter.

## 8 how to operate the digital system d-Drive®

### 8.1 common introduction

The new digital generation of piezo controllers of **piezosystem jena, d-Drive®**, combines the highest positioning accuracy with a unique handling comfort. This means that a PC can control all features and the main functions can be directly handled by the front panel. It was specially developed for multichannel functions. With the fast Digital Signal Processor (DSP) many functions are possible. We also implemented adjustable features, such as slew rate, notch filter and low pass filter. A built-in function generator offers sine, triangular, and square functions, as well as noise and sweep. The actuators are separable and exchangeable, due to integrated ID-information. The actuator type, linearity optimization, and PID-settings are stored on an ID-chip in the actuator’s plug. The calibration is done with an integrated pre-amplifier.

All these functions and the 20-bit resolution with 50 kSPS (50 kHz sampling rate) make the **d-Drive®** an excellent choice for nanopositioning and scanning applications, including scanning microscopy, metrology, surface profiling, cell tracking, nanotechnology, and many other applications.

The expandable system includes the system housing with an integrated wide range power supply in two power levels, a data exchange module with or without display, and 1 to 6 amplifier channels with digital controllers at different levels of power depending on the application.

## 8.2 system housing 84TE with power supply

### 8.2.1 function

There are two 19" 84TE system housings with two different built-in wide range power supplies with Power Factor Correction (PFC) that generate all the required high and low voltages.

On the back side of the housing there are 2 fuses, the power switch and the line filter. The short housing is prepared to operate a maximum of 6 EVD50 digital piezo amplifiers and a data exchange module with or without display. The deep housing with high power is prepared to operate a maximum of 3 EVD300 digital piezo amplifiers and a data exchange module with or without display. All modules are connected by a back plate.

### 8.2.2 technical data system housings

	84TE 315mm depth	84TE 375mm depth
Operating voltage	100 - 240 VAC / 50 - 60 Hz	
Input current	typ. 0.9A @ 115V typ. 0.5A @ 230V	typ. 1.7A @ 115V typ. 1.0A @ 230V
Power consumption	typ. 120W	typ. 230W
Primary fuses	2x T6.3A / 250V 5x20mm	
Power connector	IEC-60320 C13	
Output voltages	+145V / 300mA -30V / 300mA +15V / 500mA -15V / 500mA +5V / 2.2A	+145V / 450mA -30V / 450mA +15V / 500mA -15V / 500mA +5V / 2.2A
Dimensions (WxHxD) [mm]	450 x 150 x 360	450 x 150 x 420
Weight [kg]	4.8	7.3
Specifics	active cooling	active cooling, auto power on


table 1: technical data

### 8.3 power amplifier EVD50, EVD125, EVD300


#### 8.3.1 function

The piezo amplifiers EVDxxx are designed for low voltage piezo elements. It is integrated into a 19" housing with a width equivalent to 10 units (EVD50) or 20 units (EVD125, EVD300). It varies in output current. You can use actuators with or without an integrated measurement system (series D). It is possible to control the amplifier with an analog signal via the modulation input or digitally via CAN-Bus. The best results for resolution and accuracy can only be achieved by setting the digital target values. Because of the very low noise (lower than  $0.5\text{mV}_{\text{RMS}}$ ), the **d-Drive®** is particularly suitable for nanopositioning. The position of the actuator and all status messages are shown on the EDS2 display. If there is an EDS1 without display included, you can check the status with the command "stat *Enter*" by using a terminal program via RS232.

The new integrated digital PID-controller allows positioning of actuators which include measurement systems to eliminate drift and hysteresis. The modular construction allows a quick and easy exchange of modules and makes the system ideal for industrial purposes.

The following sketch shows the function groups of the amplifiers:

sketch 1: function groups of the amplifiers

The command for the controller consists of the analog signal applied to the MOD input, the digital value given by the OFFSET potentiometer, and the digital value from the EDS2 via CAN-Bus. The slew rate of the amplifier can be adapted to your system using the "slew rate" limit or the 4th order low pass filter, i.e. frequencies which might stimulate resonances can be suppressed before they interfere with the closed loop. The slew rate setting value is adjustable.

The PID-controller calculates the difference (err) between this input (com) and the value from the measurement system (mes):

$$\text{err} = \text{com} - \text{mes}$$

com = command  
mes = position signal

The proportional term (kp) amplifies the control deviation without reducing frequency and drives the controller until the difference between the sensed motion and the command is adjusted.

$$y_p = k_p * \text{err}$$

kp = adjustable gain  
yp = output of p-term

The P-controller can't eliminate the error completely, because it needs a deviation to drive. To minimize the permanent offset of the P-term, an integral action (I-term) is required. ki is the time constant governing the time it takes for the output to get to a certain value. For a step input it is the time taken for the output to equal the input.

$$y_i = 1/k_i * \text{err} * T_s$$

$T_s = 1 / \text{sample frequency (50 kHz)}$   
ki = time constant  
yi = output of i-term

In general the differential term combats ringing (it adds damping) and increases the reaction speed of the PI-controller:

$$y_d = k_d * 1/T_s * (\text{err} - \text{err}[n-1])$$

err[n-1] = control deviation previous sample  
kd = differential time constant  
yd = output of d-term

Be careful when using the D-term, because stochastic errors, like noise, cause extreme reactions.

The sum of these 3 terms generates a PID-control function. The output of the controller is the set value for the power amplifier to drive the piezo actuator:


$$y = y_p + y_i + y_d$$

y = output of PID-controller

In open loop the PID-controller is bridged and the command input controls the amplifier.

The following sketch describes the step response of the different terms:

jump response of loop-controllers

sketch 2: transfer function

When using a notch filter, the main resonant frequency will be damped by another order of magnitude. While using open loop you can also switch on the notch filter to reduce stimulation of the resonant frequency. Using the notch filter might increase noise.

In addition, further filters („errlpf“, „elpor“ and „tf“) are tunable available.

The filter „errlpf“ works as a low pass with selectable order to the positioning error, which is defined as the difference between the setpoint and the real position. Its order can be selected by the parameter „elpor“ with range 1st to 4th order. The positioning error represents the input value for the PID loop controller. It is preset to a default value of 180Hz with 1st order. The filter „tf“ influences the differential part of the loop control and may result in noise reduction. Typical values are  $tf=0.01 \cdot kd$ . For compatibility with earlier versions, the factory setting for  $tf$  is 0 (deactivated).

To reduce the tracking error especially during continuous operation (e.g. sinusoidal movements or similar), a pilot control (with the setpoint) can be used. The related command is „pcf,value“ with „value“ ranging from 0 (not active) to 1.0 (full setpoint feed forward). Depending on the actuator specific motion control reserve, a value of 0.75 is the default estimation. The factory setting is 0 (deactivated).

The PID-specific parameters,  $k_p$ ,  $k_i$  and  $k_d$ , are stored in the actuator's ID-chip by default. These parameters will work for most applications. Depending on your special application and load, the mechanical response can be optimized.

### 8.3.2 MOD/MON

#### modulation input: MOD

The motion of the actuator can be remotely controlled using this input. The control signal must be in the range of 0...+10 V. There is an internal addition of the MOD signal, the adjusted OFFSET potentiometer, and the present digital value. To prevent external signals from influencing any parameter, please switch off the analog modulation input when it's not used (*modon*, *x*, 0 see 8.5.3

commands). In function generator mode the modulation input is permanently switched off.

### monitor output: MON

With a special command, many different signals can be shown. The voltage range of 0...+10 V is available at this socket and can be monitored by using an oscilloscope. Pay attention to the inner resistance of the monitor output.

Command	Description	Display	Value
MON	monitor output (default = 0)	Umes/CL Udes Uact Uerr Uerrabs Uout Umes/OL	0 = position in closed loop 1 = reference input 2 = PID command voltage 3 = closed loop deviation including sign 4 = absolute closed loop deviation 5 = actuator voltage 6 = position in open loop

table 2: monitor

0 (Umes/CL) position voltage in closed loop, it means an actuator with 80 µm in closed loop generates 0...+10 V (= 8 µm/V or 0.125 V/µm)

1 (Udes) the reference input, sum of the analog signal applied to the MOD input, the digital value by the OFFSET potentiometer, and the digital value from the EDS2 via CAN-Bus

2 (Uact) output of the controller, set value for the power amplifier

3 (Uerr) closed loop control deviation including sign: the value can be positive or negative:

$$U_{err} = U_{com} - U_{mes} \quad \begin{array}{l} U_{com} = 0...+10V \\ U_{mes} = 0...+10V \\ U_{err} = -10V...+10V \end{array}$$

To generate an output of 0...+10 V the control deviation is divided by 2 and +5 V is added. To calculate the control deviation back from the measured monitor voltage, take the following formula:

$$U_{err} = (U_{mon} - 5V) * 2 \quad \begin{array}{l} U_{err} = -10V...+10V \\ U_{mon} = 0...+10V \end{array}$$

If the controller is properly adjusted the value is +5 V.

4 (Uerrabs) absolute closed loop control deviation:

$$|U_{err}| = U_{com} - U_{mes} \quad |U_{err}| = 0...+10V$$


5 (Uout) the actuator voltage -20 V...+130 V converts to 0...+10 V; to calculate the actuator voltage back from the measured monitor voltage, take the following formula:

$$U_{actuator} = (U_{mon} * 15) - 20V$$

6 (Umes/OL) The position voltage in open loop (OL) is greater than in closed loop. The monitor voltage would be smaller than 0 V and greater than +10 V. The value gets adjusted. To calculate the open loop position voltage back from the measured monitor voltage, use the following formula:

$$U_{mes}(OL) = (U_{mon} - 2.5V) * 2$$

The default setting is the position value in closed loop.

sketch 3: monitor output

### 8.3.3 pinning

#### ANALOG: SUB-D socket 9pin

pin	designation	description
1	MON+	monitor output 0 ... +10V
2	MON-	signal ground
4	MOD+	modulation input 0 ... +10V
5	MOD-	signal ground
7	TRG+	trigger output (low-active)
8	TRG-	digital ground trigger

table 4: pinning

#### PIEZO: SUB-D plug 15pin

pin	designation	description
1,2,11	AGND	analog ground
3	+15V *	operating voltage measurement +15V
4,14	GND	digital ground
5	SDA	I <sup>2</sup> C-Bus SDA
6	5Veprom *	operating voltage ID-Chip
7	Vout2	actuator voltage 2 +130 ... -20V for nanoX <sup>TM</sup>
8	Vout	actuator voltage -20 ... +130V
9	+MESS	position signal -8 ... +8V
10	-15V *	operating voltage measurement -15V
12	SCL	I <sup>2</sup> C-Bus SCL
13	DETECT	actuator detection
15	Piezomasse	actuator ground

[\*] Not for external use

table 5: pinning

The actuator is connected to this plug. The actuator voltage goes to the piezo, the position signal from the measurement system goes to the DSP, and the integrated actuator information is stored on an ID-chip which goes to the main controller. It is important that the locking screws are screwed into the front panel.

### 8.3.4 technical data EVDxxx


	EVD50	EVD125	EVD300
output power:	7.5W (max. 15W for nanoX™ actuators)	16W (max. 32W for nanoX™ actuators)	40W
output voltages:	-20V...+130V +130V...-20V for nanoX™ actuators		
output current (continuous):	50mA (2x 50mA for nanoX)	125mA (2x125mA for nanoX)	300mA (2x150mA for nanoX)
output noise:	<0.3mV <sub>RMS</sub> @ 500Hz	<0.5mV <sub>RMS</sub> @ 500Hz	
actuator / measurement connector:	15pin SUB-D plug		
monitor output:	0 ... +10V (programmable signal source)		
monitor output impedance:	50Ω		
modulation input:	0 ... +10V (disconnectable)		
modulation input impedance:	25kΩ		
trigger output (low active)	5V/0V with pull up drain resistor 240Ω		
MOD - MON connector:	SUB-D socket 9pin		
offset voltage:	-20V ... +130V via encoder		
safety functions:	short circuit proof, temperature fuse		
protection circuit:	4 fuses F50mA / 250V 5x20mm (can only be changed by authorized service personnel)		
display (LED):	green / yellow green = in operation / open loop yellow = closed loop OVL (yellow) = overload UDL (yellow) = underload		
dimensions:	3HE 10TE	3 HE 20TE	3HE 20TE
weight:	0.35kg	0.9kg	1.03kg

table 6: technical data

## 8.4 interface EDS1, EDS2

### 8.4.1 function EDS2

The EDS2 data interface enables communication between the user and the piezo drive. It displays status messages, as well as recent values of the system, on a graphic TFT-display. Specific parameters can be adjusted using the controls on the EDS2. As opposed to the EDS1, the EDS2 only offers a serial interface. It has no display or controls. All adjustments have to be done by using a PC or laptop running the terminal program. The command set is the same that is used in EDS2. For full control and comfortable adjustment we recommend using a PC or laptop running the control program. The full command set is explained in chapter titled "commands".

The POWER button switches the device on. By pressing the button for longer than 2 seconds the device gets switched off. The PC interface enables communication between the EDS and the PC or laptop. The cursor buttons navigate through the submenus. By pressing the *Enter* button one can jump into the pertaining menu. In the input settings menu, the value can be changed by using the knob. The star (\*) before the value indicates that the value was changed. By pressing the turning knob, the value is accepted. If you leave the menu without pressing the knob, the old value will be restored.

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X,Y,Z = actuator name (X=channel 0 (1. from left))  
 op/cl = open / closed loop  
 xxx.xxx = at "op" actuator voltage in Volts  
 at "cl" displacement in µm

Setup	
language	english
baud rate	115200
brightness	10
<b>back</b>	

german / english  
 9600 / 38400 / 57600 / 115200  
 0...10 brightness of TFT-display

If you place the cursor on the channel name and press Enter, the **choose menu** of the selected channel opens:

choose menu [channel]
function generator
slew rate / filter
digital controller
analog i/o
information
<b>back</b>

fct. generator [channel]	
function	sine
state	on
ampl. [%]	1.000
offset. [%]	0.000
freq. [Hz]	1.000
sym. [%]	1.000
<b>back</b>	

sine / triangle / rectangle / noise / sweep  
 on / off  
 0.0 to 100.000  
 0.0 to 100.000  
 0.1 to 9999.9  
 0.1 to 99.9

slew rate/filter [channel]	
slewr.[V/ms]	0.2
filter type	notch
state	on
freq. [Hz]	200
bandw. [Hz]	100
<b>back</b>	

0.1...100.0  
 notch / low pass  
 on / off  
 0 to 20000  
 0 to 20000

dig. controller [channel]	
KP	0.00000
KI	50.00000
KD	0.00020
<b>back</b>	

0 to 999.00000  
 0 to 999.00000  
 0 to 999.00000

analog i/o [channel]	
monitor modul.	Umess/CL on
<b>back</b>	

Umes[CL] / Usoll / Ustell / Uerr / Uerrabs / Uout / Umes[OL]  
 on / off

information [channel]	
actuator	PX200CAP
sensor	CAP
ser.no.act	64070
<b>back</b>	

actuator name  
 sensor type  
 actuator serial number

In the “information” menu all values are only read.

The “back” selection jumps back to the main menu from each sub menu.

A star [\*] before the value indicates a changed parameter. By pressing the knob the value is accepted (the star disappears). If you leave the menu without acknowledging it, the old value will be restored.

### 8.4.3 pinning EDSx

#### RS232: SUB-D socket 9pin

Pin	designation	description
2	TxD	RS232 transceiver
3	RxD	RS232 receiver
5	GND	RS232 ground

table 7: pinning

## USB 1.1: plug Series A

Pin	designation	description
1	VCC	operating voltage +5V
2	D-	data -
3	D+	data +
4	GND	operating voltage Ground

table 8: pinning

### 8.4.4 technical data EDSx

	EDS1	EDS2
Power consumption:	max. 0.9W	max. 1.2W
Voltage / current requirements:	+5V / max. 180mA	+5V / max. 230mA
Interface:	RS232 USB 1.1 CAN-Bus (internal required)	
RS232 baud rate:	115200 baud	
RS232 connector:	SUB-D socket 9pin	
CAN baud rate:	1MBaud	
Display:	<ul style="list-style-type: none"> <li>no</li> <li>Ready LED green</li> </ul>	<ul style="list-style-type: none"> <li>Full Color TFT-Display (320x240 dots) with white back light led</li> <li>Ready LED green</li> </ul>
Control elements:	<ul style="list-style-type: none"> <li>Standby/Power button</li> <li>Reset button</li> </ul>	<ul style="list-style-type: none"> <li>Standby/Power button</li> <li>Cursor buttons</li> <li>incremental encoder</li> <li>Reset button</li> </ul>
Dimensions (H W D):	3HE 24TE 160mm	
Weight:	0.2kg	0.3kg


table 9: technical data

## 8.5 communication & commands

### 8.5.1 communication via RS232

The **d-Drive** can easily communicate with your PC with the use of a terminal program (for example "hypertrm.exe" with Windows XP, please refer to your operating system's manual).

A serial connection cable (included) is required to connect to your computer.

115200 baud, 8 bit, no parity, 1 stop bit, Xon/Xoff handshake

All commands and responses of the system are encoded in ASCII-code, enabling easy reading and writing.



### 8.5.2 driver installation of USB interface

Two drivers are necessary to use the device. The first is responsible for the USB connection between PC and digital amplifier unit, the second is responsible for communication by the RS232 (COM) interface.

#### driver installation:


Please, first install the Bus driver *FTDIBUS.inf* and then the COM port driver *FTDIPOINT.inf*. The following example shows the installation procedure with Windows XP:

Please connect the PC and the d-Drive via the supplied USB cable and then turn on the system.

The device is then identified.

Please select the source.

Please select the extracted path "\psjUSB" (e.g. C:\psjUSB).



The Windows Logo testing is not relevant for this application, please continue installation.

The USB driver is now installed.


The second driver needs to be installed like the first one (COM port). It's the same procedure as before.

Please select the path "\psjUSB" (e.g. C:\psjUSB).


The Windows Logo testing is not relevant for this application, please continue installation.

The COM port driver is now installed. The installation is finished.

To use the d-Drive via the terminal program please select the COM port.


The d-Drive system must be listed in "device manager → ports".


Clicking the right button of your mouse on the device "piezojena..." enables you to change several parameters (not required).

In this example the d-Drive is connected to COM5.


**Connection setup for terminal program (Hyper Terminal):**


Please select a name for the connection (e.g. dDrive).

Please select the corresponding COM port.


The baud rate is not important. The speed is based on USB1.1.  
 Important parameters are: 8 data bits, no parity, 1 stop bit, software handshake (Xon/Xoff).

For better handling please set the following “check marks → ready”.

**Example:** operation via hyper terminal:

On top you can see the prompt “DSM Vx.xxx” after pressing *Enter*. After this, you can see the response to the “stat” command.

### 8.5.3 commands

#### global commands:

<command> *Enter*

command	description
s	shows all available commands
stat	shows the state of all amplifiers
dprpon	switch on the cyclic output of the active actuator's position value * <sup>1</sup>
dprpof	switch off the cyclic output of the active actuator's position value
dprson	switch on the automatic output of the status register when status of each module is changed
dprsof	switch off the automatic output of the status register when status of each module is changed

table 10a: commands

#### channel depended commands:

<command>,<channel>,<value> *Enter* <amplifier channel> = 0..5 (from left to right)

<command>,<channel> *Enter* without a value reads out the current value.

command	description	value
stat	status request	see table 11
mess	position value request * <sup>1</sup>	[micron]
ktemp	amplifier temperature value	[degree Celsius]
rohm	operation time of actuator since shipping	[minutes]
rgver	displays the version number of loop-controller	number
fan	switches the fan on/off	0 = off 1 = on
setf	set the output format of measurement value (mess)	0= three decimal places 1= scientific format
setg	set the output format of all floating point values except "mess"	0= three decimal places 1= scientific format
fenable	enables the actuator soft start	0= soft start disabled 1= soft start enabled
fbreak	aborts the actuator soft start	-
set	command value: actuator voltage (ol) displacement (cl)	-20 to 130.000 [V] 0 to xxx.xxx [µm] (maximum actuator displacement, see datasheet)
modon	modulation input MOD plug	0 = off 1 = on
monsrc	monitor output (0 = default)	0 = position in closed loop 1 = command value

		2 = controller output voltage 3 = closed loop deviation incl. sign 4 = absolute closed loop deviation 5 = actuator voltage 6 = position in open loop
cl	open loop / closed loop	0 = open loop 1 = closed loop
sr	slew rate	0.0000002 to 500.0 [V/ms], attributed to the modulation voltage (0 to 10V)
pcf	scale factor for pilot control	0..1 (0=deactivated, typically: 0.75)
errlpf	$f_{cut}$ of low pass filter for error value	1..10000 (typically 180) [Hz]
elpor	order of error low pass filter	1...4 (typically 1)
kp	proportional term	0 to 999.0
ki	integral term	0 to 999.0
kd	differential term	0 to 999.0
tf	filter for differential term	0..1 (default: 0; typically 0.01*kd)
notchon	notch filter	0 = off 1 = on
notchf	notch filter frequency	0 to 20000 [Hz]
notchb	bandwidth (-3dB)	0 to 20000 (max. 2 * notch_fr) [Hz]
lpon	low pass filter	0 = off 1 = on
lpf	low pass cut frequency	1 to 20000 [Hz]
sstd	set default factory settings	-----
gftk	internal function generator (see table 12)	0 = off 1 = sine 2 = triangle 3 = rectangle 4 = noise 5 = sweep
gasin	generator amplitude sine	0 to 100 [%]
gosin	amplitude offset sine	0 to 100 [%]
gfsin	generator frequency sine	0.1 to 9999.9 [Hz]
gatri	generator amplitude triangle	0 to 100 [%]
gotri	amplitude offset triangle	0 to 100 [%]
gftri	generator frequency triangle	0.1 to 9999.9 [Hz]
gstri	symmetry of triangle	0.1 to 99.9 [%] default = 50 %
garec	generator amplitude rect.	0 to 100 [%]
gorec	amplitude offset rectangle	0 to 100 [%]
gfrec	generator frequency rectangle	0.1 to 9999.9 [Hz]
gsrec	symmetry of rectangle	0.1 to 99.9 [%] default = 50 %
ganoi	generator amplitude noise	0 to 100 [%]
gonoi	amplitude offset noise	0 to 100 [%]
gaswe	generator amplitude sweep	0 to 100 [%]
goswe	amplitude offset sweep	0 to 100 [%]
gtswe	generator sweep time	0.4 to 4 [sec/decade]
sct	scan type	0 = scan function off

		1 = sine scan (one period) 2 = triangle scan (one period) 3 = sine scan (two periods) 4 = triangle scan (two periods)
ss	start scan	without value: request scan state 1 = starts scan
trgss	trigger generation stroke position start	minimum: >0.2% of total stroke to maximum: total stroke minus 0.2% of total stroke [ $\mu\text{m}$ ] or [ $\text{mrad}$ ]
trgse	trigger generation stroke position end	minimum: >0.2% of total stroke to maximum: total stroke minus 0.2% of total stroke [ $\mu\text{m}$ ] or [ $\text{mrad}$ ], always keep: $\text{trgse} > \text{trgss}$ !
trgsi	trigger generation position intervals	>0.05% of total stroke in closed loop [ $\mu\text{m}$ ] or [ $\text{mrad}$ ]
trglen	duration of trigger impulses	$n * 20\mu\text{s}$ $n=1..255$
trgedge	trigger generation edge	0= trigger generation off 1= trigger generation at rising edge 2= trigger generation falling edge 3= trigger generation at both edges
trgsrc	input source for trigger generation	0= position value (measurement value) 1= set value
trgoffs	offset for trigger values	only active by triggering on set value, in units of actuator [ $\mu\text{m}$ ] or [ $\text{mrad}$ ]
reclen	record length [samples]	0..500000
recstride	stride n for recording	1..1000 (sample interval= $n * 20\mu\text{s}$ )
recrdptr	record read pointer	0..500000
recstart	starts data recording	-
m	responds actuator position, increments record read pointer	m without parameter : m,xxxx m,1 : xxxx m,0,n : block with length n with m,xxxx m,1,n: block with length n with xxxx
u	responds actuator voltage, increments record read pointer	u without parameter : u,xxxx u,1 : xxxx u,0,n : block with length n with u,xxxx u,1,n: block with length n with xxxx
bright	brightness of TFT-display	0 – display off .. 10 – maximum display brightness

table 10b: commands

Note \*1: The modules EVDxx sends every 500ms an actual position value to the EDSx module. This value is used for display and for output via serial interface. Overuse querying this value has no effect on the refresh rate. Moving the actuator above this frequency causes alias effects, which mean the output values represents not more the real signal sequence.

### Example of command use:

An actuator with 80 $\mu\text{m}$  displacement in closed loop should have a swing in rectangle mode from

20µm to 50µm. The actuator is plugged into channel 2 (amplifier 3) --> position 3 f.l.. It should stay on 20µm for 50ms and on 50µm for 150ms. The value from measurement system is then sent to the monitor output.

commands:

- cl,2,1 *Enter*            closed loop amplifier 3 on
- gfmt,2,3 *Enter*        generator function rectangle
- gfrec,2,5 *Enter*        frequency rectangle =  $1 / (50\text{ms}+150\text{ms}) = 5\text{Hz}$
- garec,2,37.5 *Enter*    amplitude rectangle =  $(50\mu\text{m} - 20\mu\text{m}) / 80\mu\text{m} * 100\% = 37.5\%$
- gorec,2,25 *Enter*      offset rectangle =  $20\mu\text{m} / 80\mu\text{m} * 100\% = 25\%$
- gsrec,2,25 *Enter*     symmetry rectangle =  $50\text{ms} / (50\text{ms} + 150\text{ms}) * 100\% = 25\%$
- monsrc,2,0 *Enter*     set monitor source to position value in closed loop

### 8.5.4 status register

The status register is a 16bit register, in which each bit describes different properties of the amplifier or actuator. The decimal sum of all bits is the value of status register:

bit	binary	description	decimal
0	$2^0$	0 – actuator not plugged 1 – actuator plugged	0 1
2,1	$2^2, 2^1$	0, 0 – actuator without measuring system 0, 1 – strain gauge measuring system 1, 0 – capacitive measuring system 1, 1 – inductive measuring system	0 2 4 6
3	$2^3$	-----	--
4	$2^4$	0 – closed loop system 1 – open loop system	0 16
5	$2^5$	-----	--
6	$2^6$	0 – piezo voltage not enabled 1 – piezo voltage enabled	0 64
7	$2^7$	0 – open loop 1 – closed loop	0 128
8	$2^8$	-----	--
11,10,9	$2^{11}, 2^{10}, 2^9$	0, 0, 0 – generator off 0, 0, 1 – sine on 0, 1, 0 – triangle on 0, 1, 1 – rectangle on 1, 0, 0 – noise on 1, 0, 1 – sweep on	0 512 1024 1536 2048 2560
12	$2^{12}$	0 – notch filter off 1 – notch filter on	0 4096
13	$2^{13}$	0 – low pass filter off 1 – low pass filter on	0 8192
14	$2^{14}$	-----	--
15	$2^{15}$	-----	--

table 11: status register

### 8.6 function generator

Each EVDxxx amplifier in the system contains its own function generator that is able to perform multiple functions. The function generator can be programmed by using the commands detailed above. The following functions can be executed:


function type 0 = off	Sine 1	Triangle 2	Rectangle 3	Noise 4	Sweep 5
amplitude 0 to 100%	√	√	√	√	√
offset 0 to 100%	√	√	√	√	√
frequency 0.1 to 9999.9Hz	√	√	√	-	-
sweep time 0.4 to 800sec/dec	-	-	-	-	√
duty cycle 0.1 to 99.9%	-	√	√	-	-

table 12: functions and their parameters

The amplitude that has been selected is the peak to peak value. The sweep depth is fixed by 0.1Hz to 10kHz (5 decades). In the meantime, the frequency increases logarithmically. This case represents the "sweep" parameter of the increasing frequency per decade.


**Attention:**

It is not possible to calculate the waveform in real time, but it will be calculated when activate the sweep generator an stored in memory. Depending on parameter sweeptime, it can be take few seconds. When change the sweep time, initiate a new calculation of values by switching off and on the sweep generator (gfkt).


The symmetry of a signal describes, in the case of the triangular function, the ratio of rise time  $t$  to time period  $T$ . Based on this, a variation can be introduced between a fast increase and a slow decrease as well as a slow increase and a fast decrease of the function. By using the rectangular function, the duty cycle of high signal in relation to time of a period is described.

Ampl. triangle function

Ampl. rectangle function

**Attention:** Do not activate trigger function in sweep mode (trgedge =0).

### 8.7 output of trigger signals

Trigger signals allow us to get an electrical signal when a position is reached or crossed. The trigger points refer to the measured value (stroke or tilt). The trigger is low-active, i.e. a high/low edge shows that you have reached a trigger point.

The active stroke/tilt-range for trigger generation is given by "trgss" (start (lower) position) and "trgse" (end (upper) position). The interval between the trigger points is given by "trgsi". It is possible to get trigger signals at the rising edge (trgedge,1), the falling edge (trgedge,2) or both edges (trgedge,3).

To disable the trigger generation set trgedge,0. The measurement unit for trgss, trgse, trgsi is the actuator specific unit (e.g.  $\mu\text{m}$  or mrad), the length of a trigger signal can be set to multiples of 20  $\mu\text{s}$  (standard is  $1 \cdot 20 \mu\text{s}$ ).

Take care that is there no overlap between two trigger impulses. If this occurs, you must minimize the length of the impulses (trglen,1), or decrease the movement speed of the actuator.

Furthermore, you must be sure that the trigger range (between trgss and trgse) passes through plus 0.2% of the total stroke, otherwise no trigger will be generated because no change between the rising and falling edges will be detected. The points start, end and interval must be set, that the calculation  $(\text{endpoint} - \text{startpoint}) / \text{interval}$  gives a integer number  $n$  without a remainder. The counts of output trigger is given to  $n+1$ .

**Attention:** Do not activate trigger function in sweep mode.

### Example (*trgedge=1..3*):

An actuator with a closed loop stroke of 80  $\mu\text{m}$  is used as follows: rising edge (*trgedge,1*), start point 10  $\mu\text{m}$  (*trgss,10*) end point (*trgse,30*) and interval 5  $\mu\text{m}$  (*trgsi,5*), i.e. trigger points are at 10, 15, 20, 25, 30  $\mu\text{m}$ . After reaching Position 10  $\mu\text{m}$  (from a lower position; rising edge!) the trigger is set and the next trigger point will be calculated (15  $\mu\text{m}$ ), the trigger output will be reset after  $n \cdot 20 \mu\text{s}$  (*trglen,n*) and the trigger function is waiting to reach the next trigger position (15  $\mu\text{m}$ ) and so on. After reaching the 30  $\mu\text{m}$  position, the next Trigger position is set to 10  $\mu\text{m}$ , the trigger impulse will be generated when the 10  $\mu\text{m}$  position is reached from a lower position than 10  $\mu\text{m}$  (rising edge).

### Extended trigger functions (*trgedge= 4 / 5 / 7*):

The function *trgedge,4* triggers at any change of moving direction, with *trgedge,5* the output is inverted to *trgde,4*. No parameters will needet.



Abb: *trgedge,4*




Abb: *trgedge,5*

Function *trgedge,7*: at any change of moving direction a trigger with length  $n \cdot 20 \mu\text{s}$  appears (parameter  $n$ : *trglen*).

**Trigger function „moving trigger“ (trgedge,6)**

This function is parameterized by trgss, trgse, trgsi und trglen. On each period, only one trigger appears, the trigger point “moves” at every period by its increment (trgsi). After reaching the end position (trgse) the trigger point “moves” backward.

**The recognition of edges works as follows:**

During the rising edge, the highest measurement value (position) will be stored. If there is a decrease in the measurement value of more than 0.2% (of total stroke) under the stored value, a falling edge will be recognized. During the falling edge, the lowest measurement value (position) will be stored.

If there is an increase in the measurement value of more than 0.2% (of total stroke) over the stored value, a rising edge will be recognized.

This principle depends on whether or not trgss and trgse are bordered by a working area that crosses at least 0.2% (of total stroke) out of the borders.

**Triggering on setpoint:**

At very small trigger distances (near  $\leq 1/1000$  of total actuator stroke), the influence of internally electrical noise of the position value becomes noticeable. Wrong trigger outputs can be happen. A difference between the expected counts of triggers and the happened can indicate this. A solution for this problem can be triggering on setpoint. This principle uses following circumstances:

- The real positioning noise at the actuator is much less than the internally electrical noise
- The setpoint (when digitally controlled or internally generated) is noise free

The trigger generation uses the setpoint values. Consider, that the system is used in closed loop mode to reduce differences between setpoint and real position. Furthermore consider that also in closed loop mode a small difference between setpoint and real position (tracking error) is still happened. It depends on loop controller settings and moving velocity. For simple correct this error, use the parameter “trigger offset” (“trgos”). Set “trgos” to the tracking error, the the trigger outputs are shifted, so that the trigger position is well correspond to the real actuator position. The value of “trgos” is positive when the real position is following the setpoint (normally case). If the tracking error is unknown, set “trgos” to 0. On triggering by real position value, the offset is internally always set to 0.

trgss	<b>lower trigger position</b>	
-------	-------------------------------	--


	this is the lowest value where a trigger signal is generated	unit $\mu\text{m}$ or $\text{mrad}$
<i>trgse</i>	<b>upper trigger position</b> this is the highest value where a trigger signal is generated	unit $\mu\text{m}$ or $\text{mrad}$
<i>trgsi</i>	<b>trigger interval</b> this is the distance between the trigger points	unit $\mu\text{m}$ or $\text{mrad}$
<i>trgedge</i>	<b>trigger edge</b> this value depends on which direction the trigger point must be reached for trigger generation	0: trigger disabled 1: trigger at rising edge 2: trigger at falling edge 3: trigger at both edges
<i>trglen</i>	<b>trigger length</b> this values depends on the duration of a trigger impulse	$(n * 20\mu\text{s} \quad n = 1 \dots 255)$
<i>trgsrc</i>	<b>trigger source</b> selects real position value or setpoint for trigger generation	0 = real position value (measured) 1 = setpoint
<i>trgos</i>	<b>trigger offset</b> shifts the trigger points when triggering on setpoint	$[\mu\text{m}]$ or $[\text{mrad}]$

table 13: commands for trigger generation

### Hardware:

The output for trigger signals is the 9-pin D-Sub connector "ANALOG".

The pin 7 (trigger out) and pin 8 (ground) are used.

sketch 4: internal circuit

## 8.8 scan function

The scan allows the output of exactly one period of the function generator (sine or triangle). This will require that the relevant functional parameters be set (frequency, amplitude, offset or

symmetry). The function generator does not need to be activated (gfkt,0). To select a scan, set the scan type to sine (sct,1) or triangle (sct,2).


The scan starts with "ss,1". During a scan, the status can be requested with "ss<CR>". The answer "ss,2" means the scan is still running, and "ss,0" indicates the scan is complete. A running scan can be aborted by setting the status to zero (ss,0<CR>).

scan type	function	number of periods
1	sinus	1
2	triangle	1
3	sinus	2
4	triangle	2

**Example:**

scan sine , 100% amplitude, 0% offset, frequency=0,2Hz


gfsin,0.2 <CR>  
gasin,100<CR>  
gosin,0<CR>  
scan type: sct,1  
start scan: ss,1


sketch 4: output voltage of a sine scan (open loop)

**Application:**

The combined use of trigger generation and scan function permits a highly exact scan of a probe. Acceleration forces (and with it oscillations) are minimized by using a sine function, and actions can be initiated by the trigger generation at exactly defined actuator positions.

sketch 5: sine scan with trigger impulses at rising / both edges

sketch 6: triangle scan with trigger impulses at the rising edge

## 8.9 data recorder

For deeper analysis of the system behavior, a two channel data recorder can be used. Channel 1 records the position signal, channel 2 the actuator voltage. The maximum record length is 500000 samples each channel. Both channels always have the same record length, given by the command "reclen,value" [value=0..500000]. For slow or moderate movements, it is possible to reduce the maximum sample rate of 50kHz, by just recording every  $n^{\text{th}}$  value. This can be done with the command "recstride,N" [N=1..1000]. The recording duration is given by

$$t = \Delta t \cdot \text{reclen} = (20\mu\text{s} \cdot \text{recstride}) \cdot \text{reclen}.$$

There are several ways to start the data recorder:

- the "set"-command; a running recording will be aborted and a new one will started
- the start of wave generator (command "gfmt,n" , n>0)
- the start of scan function (command "ss,1")
- the command "recstart" (also without performing any movement)

The recording ends automatically after reaching the desired record length (reclen). The readout can be done independent from recording, also while recording.

Readout the data recorder:

To readout the data recorder, the commands „m“ (position value) and „u“ (actuator voltage) are used. The readout can start at any recorded position, the pointer to the start address can be set with „recrdptr,n“ [n=0..499999]. After this, the readout of the recorded position value or actuator voltage value by commands "m" or "u" can be executed. The address pointer will be automatically increased, so that a further command returns the next recorded value. This pointer increase effects both channels. Thus, to read out the same index of both channels the initial pointer value has to be set again.

There are different readout modes. When the commands "m" or "u" are sent without a parameter or with parameter=0, the answer is formatted to "command,value<CR>", when the commands are sent with parameter=1, the answer is formatted to "value<CR>". By appending a further parameter value=n, the device responds a block of n values (see examples at the end of this chapter).

### Output data format

To reduce the amount of transmitted data, a data reduced transfer with a constant data length is realized.

Counts with value range from 0 to 0xffff (16bit, hexadecimal format) will be transmitted. These counts represent a range for the position value from -30% to 130%, where 0%..100% corresponds to the closed-loop-range of the actuator and an overshoot of +/- 30% is representable.

The actuator position is calculated by:

$$\text{position} [\%] = \frac{160}{65535} \cdot \text{counts} - 30$$

For the actuator voltage range, -5%..105% is sufficient, where 0% corresponds to a voltage of -20V and 100% corresponds to +130V. Hence, the maximum representable voltage range is given to -27.5V..+137.5V.

$$\text{actuator voltage} [V] = \frac{165}{65535} \cdot \text{counts} - 75$$

A response to command m, for example, might look like: „m,b63a“. The value b63a<sub>hex</sub> represents 46650<sub>decimal</sub>, this corresponds to a position of 83,89%, based on the closed-loop-range.

$$\text{position} [\%] = \frac{160}{65535} \cdot 46650 - 30 = 83.89$$

*Example 1:*

*A jump (in open loop mode), with a record duration of t=200ms and a sample rate of Δt=100μs shall be recorded. From this follows a sample count of 2000 (you can always sample more, but minimum the desired readout length). Afterwards the position values are read out, followed by the voltage values. Thus, the minimum sample count is 2000 based on the calculation*

$$\text{rectlen}_{\min} = \frac{t}{\Delta t} = \frac{200\text{ms}}{100\mu\text{s}} = 2000$$

*Command history for the issue described above:*

```

cl,0           // open loop
set,-10       // set initial voltage to -10V

               // wait a bit, until the initial position is reached (e.g. 0.5s)

recstride,5   // 5x20μs=100μs
reclen,2000   // 200ms/100μs=2000 samples
set,110       // jump in open loop to 110V, record starts automatically

               // wait a bit, until the record is done

recrdptr,0    // rec_read_pointer (readout start address set to 0)
m             // read out 1st sample (position)
m             // read out 2nd sample (position)
...

recrdptr,0    // rec_read_pointer (readout start address set to 0)
u             // read out 1st sample (voltage)
u             // read out 2nd sample (voltage)
...
    
```

Example 2:readout modes:

u<CR>      Response:    u,value<CR>

u,0<CR>    Response:    u,value<CR>

u,1<CR>    Response:    value<CR>

u,0,3<CR>    Response:    u,value<CR> u,value<CR> u,value<CR>

u,1,4<CR>    Response:    value<CR> value<CR> value<CR> value<CR>


## 9 handling

The actuator is connected to the EVDxxx using the “PIEZO” plug. The potentiometer knob is used to adjust a static output voltage or position, depending on the operation mode (Offset). The “ANALOG” plug offers an analog interface to control the amplifier using a voltage between 0 and +10 V making arbitrary and scan functions possible. Please remove any analog input voltage from the amplifier before you turn the system on or off to avoid damage.

After power-on, the “READY” LED on the EDS lights up. When the power-on self test is finished, the display shows “main menu”. All amplifiers start their own power-on self test and initialization. When the LED's “UDL” and “OVL” stop blinking, the amplifier has initialized successfully. During this procedure all necessary parameters are transferred from the actuator's ID-chip to the amplifier and the EDS. Amplifiers without an attached actuator cannot complete this sequence. They remain blinking and cannot communicate with the EDS.

When an EDS2 is part of your system, all initialized channels are displayed. The amplifiers from left to right are shown from top to bottom, meaning that the amplifier in position1 (left) is shown in row 1 (top). For example, if positions 1, 2, and 3 are used with amplifiers, the display shows the actuator name “X” (“Y”, “Z”...) in rows 1, 2, and 3, “ol” for open loop and the actuator voltage “-20.000 V”. Rows 4, 5, and 6 remain blank, which means there was no amplifier found.

When the system starts up and no external voltage is applied to the MOD input, the actuator is in the most negative position of its deviation (see sketch). By pressing the potentiometer in, the operation mode changes to closed loop. The actuator moves to the start position of the closed loop range. Due to the different motion ranges of the actuator in open and closed loop, and because of the fast travel of the actuator to the 0  $\mu\text{m}$  position, a cracking noise can be heard. This is normal and not a malfunction! Turning the potentiometer to the right will increase the deviation until the maximum closed loop position is reached. This range will vary depending on your type of actuator.

*sketch 7: open loop / closed loop*

The yellow “OVL” LED shows an overload, while the yellow “UDL” LED shows an underload. Please try to avoid this. When an external analog voltage is applied, try to decrease or remove it. If

the LEDs still light up, check the connection between the amplifier and actuator. Additionally, a mechanical blocking or overload of the actuator can prevent actuators from reaching their position. If none of the described troubleshooting has worked, turn off the system and contact our support team. Please do not continue working after your actuator has been damaged.

**The values of calibration are valid only for a specified assembly configuration. Any change in the assembly configuration can cause the modes OVL or UDL. Please provide us with your assembly configuration in advance.**

## 10 controller adjustment

When any actuator made by **piezosystem jena** is connected to the EVDxxx amplifier, their specific values are read from the actuator's ID-chip. The Digital Signal Processor (DSP) of the amplifier is set with these values. These parameters were investigated in the **piezosystem jena** laboratory and ensure safe function of the actuator.

There are no  $k_p$ -,  $k_i$ - and  $k_d$ -values for actuators without a measurement system because the closed loop mode is impossible in this configuration. Closed loop is automatically switched off.

To adapt the controller properties to your special application, please start with the default values. First turn on the closed loop mode by pressing the OFFSET potentiometer in. Change one parameter step by step and check the result with the oscilloscope on the "ANALOG" socket of the pertaining amplifier. You must install the monitoring output to **Umes/OL** in the **analog i/o** menu.

In general, change the parameters in small steps and, depending on the actuator's reaction, slightly increase the step width. If the system begins to oscillate, switch off the closed loop immediately by pressing the OFFSET potentiometer button of the pertaining amplifier, then reset the last values entered!

First of all check the function of the notch filter. This has been factory adjusted so that the main resonant frequency is suppressed in an actuator with standard mass loading. If your application has a different mass, the resonant frequency changes (higher mass = lower frequency). The notch parameters are adjustable in the **choose menu "slewrates/filter"**. In **filter type** one has to select **notch**. The bandwidth has to be reduced.

The sweep function of the internal function generator can be used to determine the resonant frequency. Use **function generator** in the **choose menu** of the actuator. Set the amplitude to 5%, the offset to 0%, and the sweep time to 1 (1 sec/decade). Set the amplifier to open loop ("OL/CL" LED lights green). With the MOD/MON cable (MON plug) you can connect the oscilloscope to the ANALOG socket of the amplifier. The following adjustments on the oscilloscope are required: store function, time base = 0.5 sec/div, and input voltage = 0.1 V/div. The display shows the response of the actuator measured by the integrated measurement system. Before the measured curve swings out, there is a position where the amplitude is 0 Volts. This is the adjusted notch frequency. At the resonant frequency, the measured curve goes very high. At this point you have to adjust the notch frequency to the resonant frequency to linearize the curve. Using the notch filter might increase noise.

The correct  $k_i$ -value is determined as follows:

In the **function generator** menu, set rectangle to 1 Hz, amplitude to about 50%, and offset to 25%. Please set the monitoring output to **Umes/CL** in the **analog i/o** menu. The following adjustments on the oscilloscope are required: store function, time base=0.05 sec/div, and input voltage=1 V/div. Set the amplifier to closed loop ("OL/CL" LED lights yellow). The oscilloscope display shows the step response of the actuator measured by the integrated measurement system. If the slew rate is too low, please increase the  $k_i$ -value in the **digital controller** menu. The overswing after the rise time should be smaller than 1% of the total step.

In these adjustments, the actuator system can begin to oscillate in resonant frequency. Please switch off the closed loop immediately by pressing the OFFSET potentiometer button of the concerned amplifier. Then reset the last values entered! Continuous use at resonance can damage the actuator!

Now you can try to increase the slew rate in the **slewrates/filter** menu, as long as no oscillation or large overshoot occurs.

## 11 troubleshooting

Please check all cables and connections first if the system is not working properly.

error	possible solution
no reaction at switch on	Check the power cord. Check whether the main switch on the back side is on. In the main switch there are two fuses. Remove the power cord from your device and check the fuses. If one is damaged please replace, only using fuses of the same value.
UDL/OVL-LED blinks	Check the connection between actuator and amplifier. Check the actuator cable for damage.
UDL/OVL-LED light up in closed loop mode	The actuator is not able to reach the commanded position. Check the modulation input and offset. Check your actuator for mechanical blocking. Check whether your actuators move properly in open loop mode.
actuator oscillates in resonant frequency in closed loop	Check the PID-controller adjustments, reduce the ki-value, Reduce the kd-value. Check whether the function generator is still working.
read back parameters always deviate with the same level from the digital target value	Please switch off the analog modulation input when it's not used. Command: <i>modon, x, 0</i> (see 8.5.3 commands). External signals can cause interference or a command value offset.
actuator does not move in spite of the applied analog modulation signal on the MOD input	Please switch on the analog modulation input by using command <i>modon, x, 1</i> (see 8.5.3 commands).
displayed / interface values don't represent the real signal sequence	Maybe aliasing effect occurs. When the refresh rate is less than the half working frequency (under sampling), aliasing will occur. This has an effect only for display, not for the real system.

table 14: errors

## 11.1 error register

The error register is a 16bit register. Each bit describes different error. Once error has occurred the error register changes and error message as a decimal number will be issued via interface.

„?ERR,channel,error“ CR LF.

The decimal sum of all bits results the error value:

bit	description		decimal
0	I <sup>2</sup> C error	0 – no error 1 – I <sup>2</sup> C error	0 1
1	-----		--
2	temperature	0 – no error 1 – temperature out of range	0 4
3	OVL	0 – no error 1 – overload in CL	0 8
4	UDL	0 – no error 1 – underload in CL	0 16
5 – 15	-----		--

table 15: error register

For high precision positioning with nm-resolution, you must warm up the amplifier 2 hours before use. A constant temperature is beneficial. Consider that a variation of 5 Kelvin results in an elongation of 13 µm in steel with a length of 200 mm. The best results for resolution and accuracy can only be achieved by setting the digital target values. Therefore, please switch off the analog modulation input when it's not being used.

Equipment customization makes adaptations, according customer preferences, that are possible in terms of the technical threshold values e.g. the main voltage or the output voltage. Please contact our technical service department in order to find out the possibilities for your specific application. Special adaptations must be paid for by the customer.

## 12 your notes